

# Installation Guide For Exploratory Pages

---

The installation CD contains:

- a backup file (.bck) for Exploratory Pages
- DynaVox Back-up Software version 2.1
- documentation for the Exploratory Pages
- Exploratory Pages Installation Guide (this document)
- Record form

To install the Exploratory Pages you will also need:

- For your DynaVox, (English only) DynaVox Systems Software version 2.10 or greater
- For your Mac or PC computer, (English only) DynaVox System Software version 2.09 or greater, or DynaVox Backup Software version 2.1 or greater (included on the Exploratory Pages Installation CD)
- 2 MB (2,000,000) of available memory

To install the Exploratory Pages you will:

- Add a Exploratory Pages User.
- Restore the Exploratory Pages from the .bck file.
- Change into the the Exploratory Pages user.

If you do not have the software versions listed above, please call Customer Support at (800)344-1778 for directions.

## **Installation Instructions for DynaVox Backup Software version 2.1**

To install the *Exploratory Pages* user onto your DynaVox or DynaMyte you will need to have either DynaVox Systems Software version 2.09 or greater, or DynaVox Backup Software version 2.1 or greater.

If you do not have the software versions listed above, you will need to install DynaVox Backup Software version 2.1 onto your computer prior to installing Exploratory Pages. DynaVox Backup Software version 2.1 has been included on the Exploratory Pages Installation CD. It is a hybrid CD, which will work on both Mac and Windows computers.

To install DynaVox Backup Software version 2.1 use the following instructions:

### **Installing the Windows '98 or Windows '95 Version**

1. Place the CD in the computer's CD drive.
2. Select the *Start* button.
3. Highlight *Settings*.
4. Select *Control Panel*. (The Control Panel menu is displayed)
5. Double click on *Add/Remove Programs*.
6. Select the *Install* tab.
7. Select *Next*.
8. Select *Browse*.
9. Change the drive to the CD drive if not already selected.
10. Double click DynaVox Backup.
11. Double click setup.
12. Select *Finish* from the Run Installation Program.
13. Continue to click *next* until the Setup Complete Screen appears.
14. Select *Finish*. The DynaVox Backup Software 2.1 should have been successfully installed.
15. To confirm the installation and launch the software application:
  - A. Select *Start*
  - B. Select *Programs*
  - C. Select *DynaVox Backup Software*. (This will launch the software application).
16. Refer to the enclosed Installation Guide to complete the Exploratory Pages installation.
17. When you are finished installing Exploratory Pages, remove the CD from the computer and store it in a safe place.

## **Installing the Power Macintosh or iMac Version**

1. Place the CD in the computer's CD drive. An icon entitled *Exploratory Pages* will be displayed on the computer's desktop.
2. Double click on this icon to open the CD window. Search for the icon entitled "DynaVox Backup Installer."
3. Double click on the installer icon to begin the installation process.
4. Follow the prompts until the installation is complete.
5. Double click on the hard drive icon and confirm that the DynaVox Backup 2.1 icon is there.
6. Double click on the DynaVox Backup 2.1 to open the folder.
7. Double click on the DynaVox Mac Backup icon to launch the application.
8. Refer to the enclosed Installation Guide to complete the Exploratory Pages installation.
9. When you are finished installing Exploratory Pages, remove the CD from the computer and store it in a safe place.

## Installing Exploratory Pages on Your Mac or PC DynaVox System Software (English only)

Adding a new user will create the folder into which Exploratory Pages will be transferred.

### To Add a User:

1. Launch your DynaVox System Software.
2. Go to the **Setup** Menu.
3. Select **Tools**.
4. Select **User Browser**.
5. Select **Add User**. A system keyboard window appears.
6. Enter the name for the new user (e.g., Exploratory Pages).
7. Select **O.K.** The DynaVox adds a new folder containing a set of default user files. When the operation is complete, the User name appears in the viewport.
8. Go to the User Browser, use the **Up** and **Down** arrow buttons to highlight the Exploratory Pages user.
9. Select **Restore User**.
10. Select **Restore from File**.
11. Insert the Exploratory Pages installation CD.
12. Change **Current Drive** option ring to view the installation CD.
13. Locate the **explore.bck** file on the CD.
14. Select **O.K.**
15. A confirmation that Restore User function was successful will display. Press **O.K.**
16. With the **Up** and **Down** arrow buttons highlight the name of the user you just restored. Select **Change User** and press **Yes** to confirm the operation.

17. Once the user has successfully changed, the Exploratory Pages display will appear.
18. It is recommended that you do a system maintenance by pressing the **Setup** button, then pressing **System Settings**, and finally **System Maintenance**. This process may take several minutes.

## Installing Exploratory Pages on a DynaVox 3100 or DynaMyte 3100

### A. Preparing the DynaVox and DynaMyte

On your DynaVox or DynaMyte, go to the Setup menu and be sure that you are running version 2.10 software. If you are not running 2.10 software, please call a customer support representative at 1-800-344-1778 for further direction.

Next go to the Tools menu and then to the File Browser menu. The amount of available space is displayed in the upper right corner. If you have **2 MB (2,000,000)** or greater available, go to section B. If you have less than **2 MB** free on your system, complete the steps that follow to Delete a User.

**CAUTION:** If you have less than 2 MB of available memory on your system, it will be necessary to delete a user. If you have made modifications to the user that will be deleted, it is recommended that you backup the user before you delete it. For more information on backing up a user, please refer to your DynaVox Reference Manual.

The following steps **ONLY** apply if you need to increase your available memory.

To **Delete a User** from your DynaVox or DynaMyte:

1. Select **Setup**.
2. Select **Tools**.
3. Select **User Browser**.
4. Use the **Up** and **Down** arrow buttons to highlight the name of the user you want to delete.

**NOTE:** You cannot delete the active user.

5. Select **Delete User**. A confirmation window appears asking you to confirm or cancel the delete user operation.
6. Select **Yes**. The DynaVox deletes the user folder and its contents. When the operation is complete the user name is removed from the Viewport.
7. Repeat these steps until you have at least **2 MB** available on your systems (refer to the file browser).

## **B. Installing Exploratory Pages from a .bck file to your DynaVox or DynaMyte**

Adding a new user will create the folder into which Exploratory Pages will be transferred.

To **Add a User** to your DynaVox or DynaMyte:

1. Select **Setup**.
2. Select **Tools**.
3. Select **User Browser**.
4. Select **Add User**.
5. Enter the name for the new user (i.e., Exploratory Pages)
6. Select **OK**. The DynaVox adds a new user folder containing a set of default user files. When the operation is complete, the User name appears in the viewport.

The .bck files may be transferred to your DynaVox or DynaMyte using a DynaBeam or appropriate DynaVox cable in combination with any of the following software:

- DynaVox Mac (version 2.09 or above)
- DynaVox Windows (version 2.09 or above)
- DynaVox Backup Mac (2.1 or above)
- DynaVox Backup PC (2.1 or above)

**NOTE:** File transfer instructions are included in this installation guide. However, you will need to refer to your DynaVox Reference Manual to set up your computer and device for file transfer operations.

*To transfer the .bck file from your computer to the DynaVox or DynaMyte:*

7. If necessary, attach and set up the file transfer equipment (DynaBeam or transfer cables). It is recommended that equipment be shut down before attaching cables.
8. Establish the necessary file transfer setting using your DynaVox Reference Manual.
9. On your device, go to the *User Browser* menu.
10. Use the **Up** and **Down** arrow buttons to highlight the Exploratory Pages user.
11. Select **Restore User**.

12. Select **Restore via file Transfer**.

**NOTE:** A “CAUTION!” message will be presented with all transfers, explaining that this process will replace all of the existing files that have the same name as the one marked.

13. On your computer, insert the Exploratory Pages installation CD. On the CD you will find:

- 1 .bck file (Exploratory Pages – works for the DynaVox and the DynaMyte)
- 3 documents: Exploratory Pages Doc.pdf manual, Exploreinstall.pdf, Record form.pdf)
- DynaVox Backup Installer

14. Launch your DynaVox Software (2.09 or greater) or your DynaVox Backup Software (2.1 or greater).

15. Go to the *File Browser*.

**NOTE:** If you are using DynaVox Backup Software you will automatically be in the *Backup File Browser* menu.

16. Change **Current Drive** option ring to view the Exploratory Pages.

17. Locate the **explore.bck** file on the CD.

18. Access the file by pressing **Mark/Unmark Current Item**

19. Select **Transfer Marked Files**.

20. Select **Yes** to confirm operation.

**NOTE:** A “CAUTION!” message will be presented with all transfers, explaining that this process will replace all of the existing files that have the same name as the one marked.

21. Select **Start Transfer** on the computer.

22. Select **Start Transfer** on the device.

23. When the transfer is complete, press the **Done** button on both systems.

24. If necessary, disconnect the equipment from both systems.

25. With the Up and Down arrow buttons highlight the name of the user you just restored. Select **Change User** and press **Yes** to confirm the operation.
26. Once the user has successfully changed, the Exploratory Pages display will appear.
27. It is recommended that you do a system maintenance by pressing the **Setup** button, then pressing **System Settings**, and finally **System Maintenance**. This process may take several minutes.
28. Any customized changes you make from this point should be followed by backing up the user to a second source (e.g., a backup card, PC or Mac computer).